


*Patient Instructions*

Please fill all the containers.



**1** Urinate into the toilet if needed.




**2** Lift toilet seat. Place sheets of plastic wrap (e.g. Saran Wrap™) over the toilet bowl, leaving a slight dip in the centre. Place the toilet seat down.




**3** Pass stool onto the plastic wrap. Do not let urine or water touch the stool specimen.




**4** Using the spoon attached to the cap, place bloody or slimy/white (mucous) areas of the stool (if present) into each of the containers. Do not over-fill the containers.




**5** In the CLEAN container(s) (no liquid): Add stool to the "FILL LINE" and tighten cap.  
In the SAF container (with liquid): Add 2-3 spoonfuls of stool until the liquid reaches the "FILL LINE" and tighten cap.



**6** Shake the SAF container to mix the contents with the liquid.



**7** Write Name, Date of Birth and Date and Time of collection on each container.



**8** Bring the containers back to the lab immediately after collection. Delays may affect test results. If required, store the SAF and CLEAN container at room temperature. Store the CLEAN REFRIGERATE container (if one was provided) in the fridge.